

In the Claims:

Please cancel Claims 11-29, 32 and 40, without prejudice, and amend Claims 1, 3, 30, 33, 35, 37 and 38 as indicated below. The status of all claims is as follows:

1. (Currently Amended) A sizing guide for assisting in sizing a femoral prosthetic device during an orthopedic procedure, comprising:

a sizing block mountable proximate a distal end of a femur at a fixed position relative to the centerline of the femoral canal; and

a stylus member movably mounted to the sizing block, the stylus member being adjustable relative to the anterior cortex of the femur to provide a size indication for the femoral prosthetic ~~device~~; device,

wherein the stylus member includes at least one arm with a pin thereon, the pin being capable of extending towards an anterior femur, and

further wherein the pin is positioned at a fixed distance from a point of attachment between the stylus member and the sizing block.

2. (Original) The sizing guide as recited in claim 1, wherein the stylus member is rotatable with respect to the sizing block.

3. (Currently Amended) ~~The sizing guide as recited in claim 2, A~~ sizing guide for assisting in sizing a femoral prosthetic device during an orthopedic procedure, comprising:

a sizing block mountable proximate a distal end of a femur at a fixed position relative to the centerline of the femoral canal; and

a stylus member movably mounted to the sizing block, the stylus member being adjustable relative to the anterior cortex of the femur to provide a size indication for the femoral prosthetic device;

wherein the stylus member is rotatable with respect to the sizing block; and

wherein the stylus member comprises a plurality of radially extending arms having corresponding pins of differing lengths, each pin being positioned to extend toward the anterior cortex when rotated into proximity with the frame.

4. (Original) The sizing guide as recited in claim 3, wherein the position of each pin relative to the sizing block is selected to indicate a flange height of the femoral prosthetic device.

5. (Original) The sizing guide as recited in claim 1, wherein the sizing block further comprises a medial-lateral sizing feature.

6. (Original) The sizing guide as recited in claim 1, wherein the sizing block further comprises a posterior sizing feature.

7. (Original) The sizing guide as recited in claim 5, wherein the sizing block further comprises a posterior sizing feature.

8. (Original) The sizing guide as recited in claim 1, wherein the sizing block has an opening to receive a rod extending from the femoral canal.

9. (Original) The sizing guide as recited in claim 8, wherein the opening passes through the sizing block at an angle selected to generally match the angle of the femoral canal relative to the distal end of the femur.

10. (Original) The sizing guide as recited in claim 9, wherein the orientation of the opening through the sizing block is adjustable.

11-29. (Canceled)

30. (Currently Amended) A sizing guide for assisting in sizing a femoral prosthetic device during an orthopedic procedure, ~~comprising~~ comprising a sizing block mounted to a distal cut surface of a femur, at a fixed location with respect to a femoral canal of the femur, the sizing block providing an indication of a properly sized femoral prosthetic device based on the position of the sizing block relative to the femoral canal.

canal, wherein the sizing block includes an opening for slidably receiving a rod extending from the femoral canal.

31. (Original) The sizing guide as recited in claim 30, further comprising a rotatable stylus that may be adjusted to size the femoral prosthetic device based on the distance between the femoral canal and an anterior flange of the femoral prosthetic device.

32. (Canceled)

33. (Currently Amended) The sizing guide as recited in ~~claim 32,~~ claim 30, wherein the angular orientation of the opening is adjustable between at least two positions.

34. (Original) The sizing guide as recited in claim 31, further comprising a medial-lateral sizing feature to indicate medial-lateral dimensions of the femoral prosthetic device.

35. (Currently Amended) The sizing guide as recited in ~~claim 32,~~ claim 30, further comprising a posterior sizing feature to indicate posterior dimensions of the femoral prosthetic device.

36. (Original) The sizing guide as recited in claim 31, wherein the rotatable stylus comprises a plurality of arms of differing lengths.

37. (Currently Amended). ~~The sizing guide as recited in claim 36, A~~  
sizing guide for assisting in sizing a femoral prosthetic device during an orthopedic procedure, comprising:

a sizing block mounted at a fixed location with respect to a femoral canal of the femur, the sizing block providing an indication of a properly sized femoral prosthetic device based on the position of the sizing block relative to the femoral canal; and

a rotatable stylus that may be adjusted to size the femoral prosthetic device based on the distance between the femoral canal and an anterior flange of the femoral prosthetic device,

wherein the rotatable stylus includes a plurality of arms of differing lengths,  
and

wherein a pin is coupled to each arm, the pin being oriented to extend towards an anterior region of the femur when the sizing block is at the fixed location.

38. (Currently Amended) A system for sizing a prosthetic device to be located at a distal end of the femur, comprising:

means for mounting a sizing guide proximate ~~an end~~ a distal cut surface of a femur at a fixed position relative to the femoral canal; and

means for utilizing the sizing guide to perform at least one of a plurality of available femoral ~~measurements~~. measurements, wherein the means for utilizing comprises a stylus that is rotatable about an axis that passes through the means for mounting.

39. (Original) The system as recited in claim 38, wherein the means for mounting comprises a sizing block having an opening for receiving rod extending from the femoral canal.

40. (Canceled)